
tamaVet®

Canine Distemper Virus (CDV) Test (test cassette)

REF. 1172032N-05

For veterinary use only

DE	Gebrauchsanweisung	2	NL	Gebruiksaanwijzing	20
EN	Instructions for use	5	TR	Kullanma kilavuzu	23
FR	Instructions d'utilisation	8		Symbols	27
ES	Instrucciones de uso	11		Our Teams	28
IT	Istruzioni per l'uso	14			
PL	Sposób użycia	17			

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Rev. 1, 2019-03-11 RZ

1. Intended Use

The tamaVet® Canine Distemper Virus (CDV) Test is for use in the rapid detection of specific Canine Distemper Virus antigens in **conjunctiva swabs, serum, urine or faeces samples of dogs**. The test is intended for use as an aid in the diagnosis of canine distemper and is designed for professional use only.

2. Introduction and Clinical Significance

Canine distemper is one of the most significant and contagious viral diseases in dogs worldwide. The virus affects various organs at the same time, increasing the severity of the disease. *Canine distemper virus* (CDV) can infect dogs at any age. Typical distemper victims are rescued stray dogs, pet store dogs or puppies, which are usually not vaccinated or which have a questionable vaccination history. The infected dog typically infects other dogs by coughing infected respiratory secretions; but the virus can also be shed in most other bodily secretions, including urine. Canine distemper virus is therefore highly contagious.

Clinically, an infection affects the gastrointestinal tract and skin, as well as the immune and nervous systems. Transient fever and dramatic and sporadic increases in body temperature are telltale signs of infection. Generally, the symptoms of infection are highly variable and the course of the disease depends on the immune response and the viral strain. Common symptoms include various respiratory and gastrointestinal problems, but also central nervous system disorders. Canine Distemper virus causes a hardening of the footpads and nose pads, which is most evident in older dogs. Sudden death is not uncommon with this disease, and those that die as a result invariably do so from central nervous complications caused by the disease (like seizures) or from secondary bacterial infections.

Various vaccination programs against the disease have successfully controlled the virus from spreading. Vaccination of puppies begins as early as 6 weeks of age. From an epidemiological point of view it is important to note that the risk of an outbreak is minimal in areas with a high vaccination status amongst the population. Nevertheless, due to the increasing vaccination fatigue an infection cannot be fully prevented.

3. Test Principle

The tamaVet® Canine Distemper Virus (CDV) Test is a **sandwich immunoassay** which was developed for professional use by veterinarians. The assay employs specific antibodies which capture CDV antigens from the sample and thus visualise the agent. The tamaVet® Canine Distemper Virus (CDV) Test is a highly sensitive immunoassay and comprises a test strip housed in a handy test cassette.

The tamaVet® Canine Distemper Virus (CDV) Test strip consists of a sample pad, a conjugate pad, a reaction membrane and an absorption pad. The conjugate pad contains specific gold-labeled anti-CDV antibodies. The reaction membrane is precoated with additional CDV specific capture antibodies immobilised in the test line region and, as a control,

antibodies immobilised in the control line region. The absorption pad takes up the fluid applied to the test cassette.

Following the application of the sample onto the test strip, CDV antigens in the sample material react with the gold-labeled antibodies to form a complex in the conjugate pad. The fluid then migrates up the strip by capillary action. In the test line region, the formed complex is captured by the immobilised specific antibodies so that a red line appears. The presence of a red line in the test line region indicates a positive result. If the sample does not contain CDV antigens no line will form in the test line region. Additionally, a red line must form in the control line region, irrespective of the presence of CDV antigens in the specimen. The control line serves as an internal control and indicates that sufficient sample fluid has been applied to the test strip and that proper membrane wicking has occurred.

4. Reagents and Materials Supplied (5 tests/kit)

- 5 tamaVet® Canine Distemper Virus (CDV) Test cassettes
- 5 sample tubes including an applicator stick with 1mL buffer solution
- 5 pipettes (single use)
- 5 cotton swabs
- 1 package insert

5. Additional Materials Required

Timer

6. Storage & Stability

The tamaVet® Canine Distemper Virus (CDV) Tests must be stored at room temperature or refrigerated (2-30°C). The test is stable until the expiration date printed on the foil pouch. The test cassette should remain in the sealed foil pouch until use. Do not freeze tests.

7. Warnings and Precautions

- For veterinary use only.
- For professional *in-vitro* diagnostic use only.
- For single use only.
- Store tests between 2-30°C and be sure not to freeze them.
- Do not use tests beyond the expiry date.
- Do not use tests with damaged pouches.
- Perform the test within 60 minutes after opening the pouch.
- Consider the test results as invalid read after more than 10 minutes.
- Please use the stated amount of sample. An incorrect number of drops can lead to a false result.
- Place the test cassette on a smooth, horizontal surface while the test is being performed.
- Do not touch the reaction field.
- Use only the test components included in the kit. Do not mix or replace the test components with components of a different kit.
- Avoid cross contaminations by using a new sample tube, a new pipette and a new swab for each sample.
- Sample materials are potentially infectious. During the test procedure, adhere to standard guidelines for the handling

of potentially infectious material and chemical reagents. The use of protective clothing (laboratory coat, gloves, eye protection) is recommended. Materials that have come into contact with the sample material should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

The tamaVet® Canine Distemper Virus (CDV) Test was developed for the examination of **conjunctiva swabs, serum, urine and faeces samples of dogs.**

General comments:

- Perform testing immediately after specimen collection. If immediate testing is not possible, serum specimens may be stored for 3 days in a closed container at 2-8°C. For long-term storage, serum samples must be kept below -20°C. The samples should not be frozen and thawed repeatedly.
- Use only clear, non-hemolyzed serum specimens.
- Bring samples to room temperature (15-30°C) before testing. However, do not leave the samples at room temperature for prolonged periods of time.
- If specimens are to be shipped, pack them in compliance with local regulations covering the transportation of etiologic agents.
- The transportation of whole blood samples should only be performed in exceptional circumstances. The shipping of serum is preferable.
- It is important that the swab is applied with a certain pressure so that sufficient superficial conjunctiva cells are collected.
- Sampling may be unpleasant for the dog, but this is important for correct antigen detection.
- Conjunctiva specimens contain a higher CDV-antigen concentration than other sample materials in cases of an ongoing infection.
- Fresh dog faeces or pooled faeces samples are also suitable for testing.
- If testing is not to be carried out on the same day as sample collection, the faecal sample can be stored maximum for 6 days at 2°C to 8°C.
- The faecal sample must be brought to room temperature (15-30°C) before performing the test.
- Use the cotton swab to take up faecal material. It is sufficient when the tip of the swab is covered with faeces.


9. Test Procedure

1. Specimen and kit components must be brought to room temperature before performing the test.
2. Open the foil pouch and place the test cassette on an even surface.

3a. Conjunctiva swab, urine and faeces


- I. Open the sample tube and put the cotton swab with the faecal sample into the sample tube. Using the swab, stir the buffer solution in the sample tube in order to dissolve the sample in the fluid. For an optimal test result, it is recommended to wring the swab into the sample tube by squeezing the tube with the swab in the middle several times to dissolve the sample.
- II. Remove the swab and close the sample tube. Shake the sample tube to prompt the buffer solution to react with the sample. Break off the the tip on the lid of the sample tube.
- III. Dispense **3 drops of the sample solution** into the sample well. Hold the sample tube vertically with the lid down over the sample well and squeeze the tube gently. Start the timer.

If the liquid is not migrating sufficiently along the strip after 60 seconds, add an additional drop of the sample solution into the sample well and/or push with the tip of the pipette into the sample well to reactivate the flow.

3b. Serum

- I. Holding the pipette vertically, dispense **1 drop (30µL) of serum** into the sample well. Wait a couple of seconds until the sample fluid is completely absorbed.
- II. Open the sample tube by breaking off the tip on the lid and add **2 drops of the buffer solution**, which is in the sample tube, into the sample well. Hold the buffer tube vertically with the lid down over the sample well and squeeze the tube gently. Start the timer.

Avoid the formation of air bubbles. If the liquid is not migrating sufficiently along the strip after 60 seconds, add an additional drop of the buffer solution into the sample well and/or push gently with the upper tip of the swab into the sample well to reactivate the flow.

4. Read the results 10 minutes after starting the timer. **Do not interpret the result after more than 10 minutes.** The test and the result are considered as invalid after this period.

10. Result Interpretation

Positive result

The test is positive when the control line (C) and the test line (T) are visible in the reaction area. CDV antigens have been detected in the sample material. If a weakly defined test line appears, the test result is nevertheless positive.

Negative result

The test is negative when only the control line (C) is visible in the reaction area. No CDV antigens have been detected in the sample material.

Invalid result

If no control line is visible after the test was performed, the test is invalid. All results, that show no control line after the denoted time, must be discarded. Review the test procedure and repeat the test with a new test cassette.

If problems continue, discontinue using the test kit and contact your distributor.

11. Quality Control

An internal procedural control is included in the test cassette:

The formation of a red line in the control line region (C) is considered as internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Limitations

- Test results must always be interpreted in the light of all the medical findings and never individually. All test results which do not correlate with the clinical picture must be repeated. Positive results detected with the tamaVet® veterinary rapid tests should be confirmed using another analytical method.
- Due to the CDV infection cycle, the virus can usually be detected from the beginning of the viremia phase (presence of infectious virus particles in blood), approximately 7 days after the initial infection.
- It is recommended to wait 7 days after a vaccination before testing, since a cross reaction cannot fully be ruled out.
- There is a possibility that the test result is distorted by technical errors, errors in the experimental procedure or by substances or factors that influence the test and are not mentioned here.

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nal von minden GmbH