

MOBILE LIVING MADE EASY.

DOMETIC DS 301H

NOISELESS MEDICAL REFRIGERATOR, RIGHT HINGED, 27 L



**MEDICATION
COOLING**
HEALTH CARE
SUITABLE

**TEMPERATURE
CONTROL**
DIGITAL
THERMOMETER

**SAFE
STORAGE**
AUDIBLE
WARNING

SAFE AND SECURE

The perfect solution for hospitals and health centres to keep heat sensitive medicines between +2 °C and +8 °C. Integrated fan for constant temperature and audible warning if door left open for 20 s.

With its reliably constant temperature courtesy of an integrated fan, the Dometic DS 301H is the ideal solution for hospitals and health centres wishing to store heat sensitive medicines. Contents are cooled to between +2 °C and +8 °C. An audible warning will be triggered if the door is left open for more than 20 s. The miniCool DS 301H has a door lock for secure storage and drawers with flexible compartments for organised arrangement of medicines.

- Keeps heat sensitive medicine between +2 °C and +8 °C via noiseless absorption cooling.
- An exterior thermometer displays the inside cooling performance. Integrated fan for reliably constant temperature
- Features an audible warning system if the door is open for more than 20 s?
- Storage volume 27 l
- The door has a lock for secure storage

SPECIFICATIONS

GENERAL

SKU number 9105203202
Model DS301HFS
Model on label DS301H
Scope of delivery 1 fridge
Refrigeration technology Absorption

DIMENSIONS

Dimensions product depth 393 mm
Dimensions product height 580 mm
Dimensions product width 422 mm
Gross volume (DIN15502) 28 l
Net weight 15 kg
 27 l

ELECTRICAL

Input frequency 50/60 Hz
Input voltage (AC) 230 V
Rated input power (AC) 65 W

PERFORMANCE

Average energy consumption at AC (DIN15502)
 0,86 kWh/24h
Climate class (DIN15502) SN
Noise (DIN60704-1 & 60704-2-14) 0 dB(A)
Refrigerant, type NH3
 101 g

ADDITIONAL

Colour White
Certificates CB,CE
Hinge type Right hinged
Interior light Yes, LED
Thermostat control Electronic Thermostat

LOGISTICS

Dimensions parcel depth 436 mm
Dimensions parcel height 630 mm
Dimensions parcel width 465 mm
Parcel weight 16 kg
EAN-13 5999024847737

Mobile living made easy.

