

## Product Specification Sheet – eks Instrument Brushes



### Article numbers:

- ek 11-80 – Brush with stainless steel bristles
- ek 11-82 – Brush with brass bristles
- ek 11-84 – Brush with nylon bristles
- ek 11-86 – Brush with nylon bristles, double sided

### Characteristics

- Ideal for cleaning of all types of instruments and burs
- Instrument Brushes with nylon bristles are suitable for autoclave up to 134°C
- Ergonomic handle

### Instructions for use

- Please test the eks Instrument Brushes before use at a non-visible area. We cannot accept any liability for any damages caused as a result.

### Technical Specifications

- Handle material : plastic
- Color: blue
- Bristles materials : stainless steel, brass, nylon
- Length: approx. 18 cm

The content of this manual is based on our current knowledge and experience in the development and manufacture of lubricants. Because of the complexity of tribological systems, the effect of our products depends on many parameters, which we cannot assess and which influence we cannot evaluate. For this reason general statements about the function of our products are not possible. The information in this manual, therefore, contains non-binding guidelines, which should give the technical trained reader information on possible applications. The information in this manual does not include property assurances or warranties or guarantees to the properties or suitability of this product in a specific application. Prior to its use it is absolutely necessary to test this product in the application to ensure that the product and its use is safe, economical and fully suitable. It should proceed with due diligence. We reserve the right to change the information in this manual at any time and without notice. All previous versions of this manual are no longer valid. Updating : 2. June 2016

