

SURGICAL GOWN “REINFORCED” L



DESCRIPTION

The robust **Teqler surgical gown “Reinforced”** is individually sterile packed and available in three different sizes. It is particularly suitable for surgical procedures with high liquid exposure

PRODUCT DETAILS

- Surgical gown with reinforcements in forearm and chest areas
- Size. L
- Liquid repellent - ideally suited for operations with high liquid exposure
- High resistance to seepage and germ transmigration
- Material: SMMS (Spunbond - Meltblown - Meltblown - Spunbond)
- Knitted cuffs on the arms
- Velcro closure on the neck
- Overall liquid resistant and liquid-proof in the reinforced areas
- With straps for individual adjustments
- Abrasion resistant, low linting
- 100% fibreglass free
- Colour: green
- Individually sterile packed
- 1 surgical gown, delivered with 2 sterile, non-woven hand towels

Overview

Article Number:
T135281L



Category:
Protective clothing OR

EAN:
4260306776245



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

