

## REAGENT BOTTLE | 250 ML



### DESCRIPTION

Teqler reagent bottle made of temperature and chemical resistant 3.3 borosilicate glass. Suited for storage as well as for setting and mixing solutions in the lab or practice. This glass reagent bottle is delivered with a pouring edge and a tightly sealing screw cap.

### PRODUCT DETAILS

- Glass reagent bottle
- Made of 3.3 borosilicate glass
- With pouring edge and screw cap
- Chemical and temperature tolerant (up to 140°C)
- With DIN thread GL 45
- Easily readable scaling
- Various sizes available
- 1 piece

### Overview

Article Number:  
**T135849**



Category:  
**Lab Equipment**

EAN:  
**4260547369282**



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  [info@netmed.lu](mailto:info@netmed.lu)

Social Media

