

## PLASTIC FUNNEL



### DESCRIPTION

Teqler plastic funnel for filling vessels and bottles. These funnels have an eyelet for hanging and are available in various sizes.

### PRODUCT DETAILS

- Teqler plastic funnel
- Made of transparent polypropylene
- With eyelet for hanging
- With size marking
- Various sizes available
- 1 piece

#### Fill Opening Spout Length Spout Opening

50 mm	40 mm	6 mm
75 mm	60 mm	10 mm
100 mm	80 mm	14 mm

### Overview




Article Number:  
**T135863-T135865**

Category:  
**Lab Equipment**



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  [info@netmed.lu](mailto:info@netmed.lu)

Social Media

