

## SPECIFICATION

### SERVOspin Smart Micro Hematocrit Centrifuge



Compact micro hematocrit centrifuge with interchangeable rotors for micro vessels.

The SERVOspin SMARTpro Micro Haematocrit Centrifuge is a high speed table centrifuge for the determination of the hematocrit value. Up to now the determination of the hematocrit value (percentage by volume of cells in blood) by centrifugation as gold standard (DIN 58933-1). Equipped with a low-noise motor, the SERVOspin SMARTpro Micro Haematocrit Centrifuge accelerates the rotor to a speed of 12.000 rpm which corresponds to a maximum RCF of 9660 G. The clear design allows intuitive handling by the user. The locking of the lid is monitored electronically and is coupled with an automatic emergency brake function. This offers maximum comfort and safety for the user. A 12-times hematocrit rotor with HKT evaluation card, a dispenser tube with 100 capillaries, a plate with sealing wax, two exchangeable rotors for mini centrifuge tubes, as well as a detailed operating manual are included in the scope of delivery.

#### All functions at a glance

- Maximum speed 12.000 min<sup>-1</sup> (RCF 9660 G)
- Simple and intuitive operation
- Running time and speed adjustable
- Low noise engine
- Electronically monitored lid locking system
- 24 months warranty
- 

Technical data:

Speed: 500 - 12000 rpm

Speed steps: 500 rpm

Max. RCF = centrifugal force: 9660 G

Running time: 1 to 99 min

Rotor equipment: HKT rotor, angle rotor, strip rotor, adapter set

Dimensions: 264 x 229 x 120 mm

Mains voltage: 85-264 Volt/ 50/60 Hz/ 45 Watt

Weight: 2 kg

WEEE Registrierungs-Nr. des Herstellers: DE74434443

## SPECIFICATION

### SERVOspin Smart Micro Hematocrit Centrifuge

Type	PZN	REF	Unit
SERVOspin SMARTpro Mikro-Hämatokrit Centrifuge		H7 680	1 piece

Product	REF	Quantity	Unit price	Total

\_\_\_\_\_ plus tax

Personal details:

Customer number: \_\_\_\_\_

Company: \_\_\_\_\_

First name: Last name: \_\_\_\_\_

E-mail: \_\_\_\_\_

Street, street#: \_\_\_\_\_

Postal code, city: \_\_\_\_\_

I hereby place a binding order:

Company stamp / Address

\_\_\_\_\_ Date / Sign

You are ordering under our terms and conditions, our valid terms and conditions can be found at <https://www.servoprax.de/en/terms-and-conditions/>