

HEAD FIXTURE FOR SPINEBOARDS



DESCRIPTION

With this head fixation a complete immobilization of the head can be achieved in injured patients. The supplied straps allow secure attachment to the Teqler Spineboard.

PRODUCT DETAILS

- Head fixation for Spineboards
- Adjustable head strap for secure fixation
- Waterproof, easy-to-clean plastic
- Robust and durable
- Fastening straps included in scope of delivery
- Size: 40 x 25,5cm
- 1 piece

Overview



Article Number:
T138245

Category:
Emergency Equipment

EAN:
4260306772025

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

