

TRANSFER NEEDLE



DESCRIPTION

The Teqler transfer needle consists of 2 sharp puncture spikes and a grip plate that facilitates handling. The transfer set can be used, for example, to mix solutions and medicines or to dissolve lyophilisates.

PRODUCT FEATURES:

- Transfer cannula
- Sharp, slim puncture mandrels
- With grip plate in the middle
- Total length: 6.3 cm
- Individually sterile packed
- Scope of delivery: 1 piece

Overview

Article Number:
T138323



Category:
**Infusion | Injection |
Syringes**

EAN:
4260626263838



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

