


DIS[INJECT

ROTATABLE HEAD POLES



injection




surgery



orthopedics



hospitalization



capture



protection



hygiene



livestock



grooming


génia
the veterinarian's preference

CAPTURE POLES

■ ROTATABLE HEAD POLE

Pole with a rotatable head, professional quality, for capturing dogs of all sizes.

- 2 non-slip grip handles **A**
- Aircraft-quality 22 mm diameter, 1.7 mm thick aluminium reinforced tube with protective coating **C**
- 3 mm multi-strand stainless steel cable with plastic coating. Capacity of up to 450 kg.
- Ball handle for locking / unlocking the cable (widening of the capture loop through an automatic cable return) **B**
- Anti-strangulation rotatable head **D**




Code	Model	Length	Weight
221430	Rotatable head pole, reinforced tube	120 cm	0.96 kg
221432	Rotatable head pole, reinforced tube	150 cm	1.125 kg

■ TELESCOPIC POLE

Professional quality telescopic pole with rotatable head for capturing dogs of all sizes.

- Identical to rotatable head pole, with an extension capacity from 120 to 180 cm.
- Locking ring **E**



Code	Model	Length	Weight
221434	Telescopic pole	120-180 cm	1.16 kg



Repair on request.



CAPTURE