

# InSight V5<sup>®</sup>

## Veterinary Haematology Analyser



**5 PART DIFFERENTIAL  
WITH LASER TECHNOLOGY.**



**PRAXISDIENST**  
Medical Supplies since 1953  
[Order here!](#)

**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

# InSight V5<sup>®</sup>

**A fully automated multi-species,  
5 Part Differential Haematology  
Analyser with Laser Technology**

## FEATURES

- Quick and easy to use
- Small sample volume - 15µl whole blood
- 60 second analysis time
- 5-part differential, 23 parameters
- Species specific reference ranges
- Independent channel for Basophil measurement
- Powerful capability to flag abnormal cells
- Fully automated
- Touch screen
- Low maintenance
- Results transmitted via InSight p-LINK to Practice Management Systems



## SPECIFICATIONS

### Principles

- Tri-angle laser scatter method for classifying and counting WBC
- Impedance method for RBC and PLT counting
- Cyanide free reagent and colourimetric method for haemoglobin concentration testing

### Parameters

WBC	LYM#	HGB	HCT
NEU%	EOS%	MCV	PLT
NEU#	EOS#	MCH	MPV
MON%	BASO%	MCHC	PDW
MON#	BASO#	RDW-CV	PCT
LYM%	RBC	RDW-SD	

### Species :

Dog, Cat, Rabbit, Horse, Cow, Guinea Pig, Mouse, Pig, Rat & Sheep

### User Interface

- Easy to use, menu driven user interface
- Screen Type - LED backlit, colour graphics LCD

### Connection Interfaces

- One LAN interface
- Printer Options - Built-in thermal printer or connect to external inkjet and laser printers
- 4 x USB interfaces

### Graphic Parameters

- Four 2D scattergrams
- Three histograms

**Integrates with  
Practice Management Systems**



### Software Languages

English, Chinese, Spanish, French, Russian, Portuguese

### Dimensions

- Width 36cm (14.3")
- Depth 43cm (17")
- Height 50cm (19.6")
- Weight - 28kg (61.7lbs)
- Power Supply - 100-240V AC

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**