

## SHARP CURETTE 16,5 CM



### DESCRIPTION

The **sharp curette from Teqler** is a double-headed instrument equipped with 2 functional heads of different sizes. The sharp curette, which is also known as a Partsch bone curette, is used for the removal of tissue overgrowth or necrotic tissue, for example. The curette is made of matted steel and is clearly labelled as a single-use instrument to prevent any confusion with reusable instruments.

### PRODUCT DETAILS

- Teqler sharp curette
- Length: 16.5 cm
- Material: steel, matted
- Double-headed instrument with 2 functional heads
- Colour-Coded
- With labels
- Fixed on cardboard dispenser sheet
- 3 year sterile shelf-life
- Individually packed sterile
- Dispenser box with 20 single-use instruments

### Overview

Article Number:  
**T710114**



Category:  
**Clamps | Holders | Pliers**

EAN:  
**4260306770144**



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  [info@netmed.lu](mailto:info@netmed.lu)

Social Media

