

# Silamat<sup>®</sup> S6

Form and function –  
perfectly mixed



- Aesthetic design
- Homogeneous mixing results
- Universal use

# Silamat® S6

## Form and function – perfectly mixed

The Silamat S6 is a modern universal mixing device for amalgam, glass ionomer cements and other predosed dental materials in capsules and injection capsules.



### Form and function

The aesthetic design is not only attractive but also functional. Thanks to its rounded shape and the removable cover the unit can be cleaned easily.



### Homogeneous mixing results

The patented and time-tested Silamat® technology with the figure-eight mixing pattern and above-average rotation speed of 4500 rpm help to achieve reliable and homogeneous mixing results.



### Universal use

The Silamat S6 can be used to mix the contents of all popular capsules and injection capsules in various shapes and sizes.



### Easy handling

Simply load the capsule from above with one swift hand movement, set the mixing time, start, and you're done!

## Silamat S6

## Art. No. 602 286

### TECHNICAL DATA

Rotations per minute  
Operating voltage  
Dimensions  
Warranty

4500 rpm  
100-240 VAC  
LxWxH = 240x230x170  
3 years

Descriptions and data constitute no warranty of attributes.  
Printed in Switzerland  
© Ivoclar Vivadent AG, Schaan / Liechtenstein  
611508/0208/e/RDV

Manufacturer and distributor  
**Ivoclar Vivadent AG**  
Bendererstr. 2  
FL-9494 Schaan  
Principality of Liechtenstein  
Tel +423 / 235 35 35  
Fax +423 / 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

  
passion vision innovation